

## Routine Equipment Change List

Transducers	96hrs (change at cannula – not at the tap, make sure all admin sets and extensions are changed)
IV administration sets for infusions that are continuously used(except 10% or 20% Glucose)	96hrs
<b>IV octopus</b>	<b>TBC - 96 hours</b>
Lines that are intermittent (single use)	Change after infusion ends
<b>10% or 20% Glucose,</b>	<b>24hrs</b>
Propofol administration sets	12hrs (6hrs if infusion is via syringe driver) 12 hrs
IV drug mixing spiked bag	24hrs
Blood administration sets	After completion of max 2 units of transfusion
Enteral Feeding sets	24hrs
Urinary Catheter bags	7 days
HME	24 hrs
Peripheral IV cannulae	Remove if signs of inflammation or VIP>2 Label dressing with date of line insertion.
CVC dressings.	7 days (or earlier if wet, soiled or loose) Date with dressing change.
Closed suction catheters	72hrs
Suction tubing & yankauer	24hrs
LiDCO™ sensors/park + ride.	48hrs
cRRT circuit	72hrs
Nebuliser sets	24hrs each morning Wash with sterile water & air dry
Sterile water bottles	24hrs

**On admission- MRSA SCREENING: Nasal (both nostrils, 1 swab), Groin (both, 1 swab)  
Pseudomonas SCREENING (Throat, Perineum, Sputum, CSU, Wounds, but NOT nose or axilla)  
Weekly MRSA SCREENING on ALL patients. Send Monday am (Tuesday if bank holiday weekend).**

**Please label all lines with date and time**